

## On-Line Panel Mounted pH Indicator MS pH 19



## **INTRODUCTION**

MicroSet Digital pH meter Type MS pH 19 is a precise and compact pH Indicating System using latest technology of electronics for direct measurement of pH. The high input impedance, accuracy and the stability of its measuring circuit enable it to make the measurements with precision and reproducibility.

## **SPECIFICATION**

Range	: 0~14 pH
Accuracy	: ±0.1 pH
Display	: 31/2 Digit, 14 mm Red LED
Repeatability	: ± 1 LSB
Power Supply	: 230 V AC, 50-60 Hz
Enclosure	: ABS
Size	: 96 x 96 mm

Supplied with pH buffer solution of 4.0 & 7.0 1 no each